

The material in this folder is divided into three sections. The first section provides basic information on the biochemistry of DNA and is also intended to acquaint the students with the principles and techniques of electrophoresis. The second part shows students how to apply what they have learned to perform one or more exercises in modern biology. Each exercise consists of a background information section, an experimental procedure and study questions. Before the laboratory, the student should read the background information section, study the directions for doing the experiment and understand the reason for each step in the procedure. The instructor guide for experiments 1001 - 1004 is attached to the front of the student's instruction sheets. This guide contains a listing of the chemicals provided with each experiment, procedures for preparing solutions, typical experimental results, and answers to study questions.

John N. Anderson, Ph.D.
 Professor of Biology
 Department of Biological Sciences
 Purdue University
 West Lafayette, Indiana 47907

© 1993 by John N. Anderson, All rights reserved.

NOTES:

1. Permission is granted to reproduce the written material in this folder one time for educational purposes only.
2. Electrophoresis Package 3/4 is required to perform the exercises described in this folder.




Table of Contents

Part A. Background Information

I. Nucleic Acid Structure and Function: A Review of the Basics	
Nucleotides-Building Blocks of Nucleic Acids.....	1
The Polynucleotide Chain.....	3
DNA Structure.....	4
Biological Role of DNA and RNA.....	5
Analysis of Specific Genes and Recombinant DNA Technology.....	7
II. General Description of Agarose Gel Electrophoresis	
III. Practical Aspects of Electrophoresis	
Electrophoresis Equipment.....	12
Electrophoresis Chemicals.....	14
IV. Electrophoresis Procedures	
Pouring the Agarose Gels.....	15
Sample Application.....	16
Electrophoresis.....	17
Staining.....	18

Part A. Background Information

I. Nucleic Acids: A Review of the Basics

The concept that chromosomal units known as genes transmit heritable information from parent to offspring was founded in the late 19th century. However, a description of genes in terms of their unique structural and functional properties is relatively new. We now know that genes are composed of a type of nucleic acid called deoxyribonucleic acid (DNA). The DNA molecule not only directs its own reproduction but also stores all the information that determines the types of proteins produced during the lifetime of an organism. In so doing, DNA orchestrates the complex reactions and structures characteristic of an organism and its offspring. Ribonucleic acid (RNA), the second major category of nucleic acids, is involved principally in the transmission of genetic information and in protein production. The structure and function of DNA and RNA can most easily be understood by examining the chemical composition of the nucleic acids.


Nucleotides – Building Blocks of Nucleic Acids

Under the proper conditions, nucleic acids can be broken down to low-molecular-weight products of three types: a pentose (or 5 carbon) sugar, purines and pyrimidines; and phosphoric acid (Figure 1). The phosphate group is responsible for the strong negative charge of nucleic acids. The pentose sugar from RNA is always ribose and that from DNA is 2-deoxyribose. These sugars differ only by the presence or absence of a hydroxyl group on carbon 2 (so-called 2'). The numbers assigned to the five carbon atoms are shown in Figure 1. The purines and pyrimidines are often called nitrogenous bases (or, simply, bases). The major purine bases in DNA and RNA are adenine (A) and guanine (G), and the major pyrimidines in DNA are cytosine (C) and thymine (T). RNA contains the base uracil (U) in place of thymine. The sugars and phosphates are readily soluble in water. That is, they are hydrophilic. In contrast, the bases are hydrophobic in that they display limited solubility in water. As will be discussed below, these differences in water solubility are extremely important for the structure of the DNA molecule.

A nucleotide consists of a pentose sugar, a nitrogenous base and a phosphate group structured as shown below. The high-energy storage compound, adenosine triphosphate (ATP), is a well known nucleotide found in biological systems.




Figure 1. The Nucleotide Components




The Polynucleotide Chain

Nucleic acids are polynucleotides and have the general structure shown below:



A polynucleotide is composed of repeating nucleotide units linked into chains by phosphodiester bonds that join the 5' carbon of one ribose or deoxyribose group to the 3' carbon of the next sugar (Figure 2). The sequence or order of nucleotides in a polynucleotide chain is often abbreviated by a 1-letter code (e.g., G-C-A-T-A) with the 5' end of the chain written at the left. A typical RNA molecule is a single-stranded polynucleotide chain. As will be described below, DNA usually contains two polynucleotide strands coiled around one another to form a double-stranded helix. The number of nucleotide units in a nucleic acid chain varies tremendously depending on the nucleic acid type. For example, each chromosome from a higher organism is thought to contain a single, very long DNA molecule. A DNA molecule from the largest human chromosome is composed of approximately 5.4×10^8 nucleotides, which corresponds to a molecular weight of the order of 10^{11} and a length of about 4 cm. On the other hand, transfer RNA molecules generally contain only 70-80 nucleotides.

Figure 2. Structure of the Polynucleotide Chain




3

DNA Structure

In early physical studies of DNA, a variety of experiments indicated that DNA molecules occur in long helices with each helix being formed from two or more polynucleotide chains bound side by side. Chemical analyses also demonstrated that the phosphate groups were on the outside of the helix and that the number of A and T residues in DNA were always equal, as were those of G and C. With these facts in mind, Watson and Crick in 1953 proposed that the DNA molecule actually consists of two polynucleotide chains coiled around the same axis to form a double helix (Figure 3). In this model, the hydrophilic sugar-phosphate groups follow the outer edges of the molecule where they can interact with water. The hydrophobic bases face inward toward each other in the molecule's center and thus avoid contact with water. The two polynucleotide strands run in opposite directions (they are antiparallel) and are held together primarily by hydrogen and hydrophobic bonding between the bases, where A is always paired with T, and G with C. These complementary bases have an affinity for each other such that, when they are paired, they contribute to the overall stability of the DNA helix. Because of this complementary base-pairing, the sequence of bases in one polynucleotide chain determines the sequence in the other. For example, if the bases along one strand are arranged in the order T-G-C-T-A-G, the opposite bases on the complementary strand will be A-C-G-A-T-C. This fact is of extreme biological significance because it explains how a DNA helix in the chromosome directs the formation of copies of itself and directs the formation of RNA molecules with its specific informational content. The B-form DNA shown in Figure 3 is the most common of the DNA types. It is a right-handed double helix and contains about 10.5 nucleotide pairs per helical turn.

Figure 3. Double-Helical Structure of Common B-DNA



4


Biological Role of DNA and RNA

DNA is an information molecule with two general functions (Figure 4). First, DNA plays a central role in the propagation of the species and the determination of the heritable characteristics of the cell and its descendants. Prior to the time of each cell division, the two strands of the DNA helix separate from one another and each serves as a pattern or template for the synthesis of a new, complementary chain. This process of DNA biosynthesis is known as replication. One of the double helices formed is then transmitted to one daughter cell, and one to the other. Although the principle underlying DNA replication is straightforward, the actual mechanism responsible for the replication process in the cell involves an array of enzymes and regulatory proteins.

The informational content of DNA also determines the types of proteins that are produced by a cell. In this manner, the DNA molecule functions as a blueprint for all cellular processes that go on during the lifetime of an organism. In the first step along the information pathway from DNA to protein, a segment of DNA is copied into a complementary strand of messenger RNA (mRNA) by a process known as transcription. Transcription begins when an enzyme called RNA polymerase binds to a specific sequence on the DNA known as the promoter. At this site, the enzyme unwinds a small segment of double helix, exposing the bases of the two single strands of the DNA molecule. One of these strands is then transcribed. As the polymerase travels along the DNA, ribonucleotides with bases complementary to the DNA are added to the growing chain. For example, the DNA segment C-G-T-A-T-G is transcribed into G-C-A-U-A-C in the mRNA. Each sequence of three nucleotides in the mRNA is called a "codon," and codes for one amino acid. Since most polypeptide chains contain between 100 to 1,000 amino acids, an mRNA must be at least 300 to 3000 nucleotides long. Therefore, a gene that codes for a polypeptide chain must contain at least 300 to 3,000 base pairs.

The translation of mRNA into protein is a complex process that occurs on particles called ribosomes. This process requires ribosomal RNAs (rRNAs) and transfer RNAs (tRNAs). These RNA species do not specify proteins themselves but rather take part in decoding the information carried by the mRNAs. At one end of each tRNA molecule is a nucleotide triplet called the "anticodon," which is complementary to an mRNA codon. A specific amino acid is bound to the opposite end of each tRNA molecule. On the ribosome, tRNAs carrying amino acids associate with the mRNA by way of complementary base-pairing at the anticodon-codon sequences. As the ribosome moves along the mRNA, the amino acids carried by the tRNAs are linked to the growing polypeptide chain. In this manner, the order of codons along the mRNA directs the amino acid sequence of a polypeptide chain. The translation process occurs on the surface of the ribosomes. These particles, which are composed of rRNA and ribosomal proteins, serve to bring together the mRNAs, the tRNA and other factors that are required for protein production.

Figure 4. Molecular Information Transfer: Complementary Base Pairing



Analysis of Specific Genes and Recombinant DNA Technology


A key to one of life's great mysteries was discovered in 1953 when the double-helical structure of DNA was perceived by Watson and Crick. Elucidation of the basic mechanisms of replication, transcription and translation quickly followed, and by the early 1960's, the model shown in Figure 4 was generally accepted by most biologists. However, genes from higher organisms resisted detailed analysis until the mid 1970's because of the complexity of the DNA in eukaryotic organisms; a vertebrate cell contains enough DNA to code for more than 100,000 proteins. In order to study the structure and function of a single protein coding gene, the gene must be prepared in a purified form. The isolation of a specific gene from cellular DNA by conventional biochemical procedures is not practical because of the magnitude of the purification required (usually 100,000-fold) and because the procedures would necessitate the use of a large quantity of starting cellular DNA. Herein lies the major use of recombinant DNA technology, for it permits the amplification and isolation of specific genes by relatively simple procedures. A basic understanding of these procedures requires a description of an interesting feature of bacterial physiology.

Plasmids are small circular DNA molecules that exist apart from the chromosomes in most bacterial species. Under normal circumstances, plasmids are not essential for survival of the host bacteria. However, many plasmids contain genes that enable the bacteria to survive and prosper in certain environments. For example, some plasmids carry one or more genes that confer resistance to antibiotics. A bacterial cell containing such a plasmid can live and multiply in the presence of the drug. Indeed, antibiotic-resistant *E. coli* isolated in many parts of the world contain plasmids that carry the genetic information for protein products that interfere with the action of many different antibiotics.

In the laboratory, plasmids can be introduced into living bacterial cells by a process known as transformation. When bacteria are placed in a solution of calcium chloride, they acquire the ability to take in plasmid DNA molecules. As illustrated in Figure 5, this procedure provides a means for preparing large amounts of specific plasmid DNA since one transformed cell gives rise to a clone of cells that contains exact replicas of the parent plasmid DNA molecule. Following growth of the bacteria in the presence of the antibiotic, the plasmid DNA can readily be isolated from the bacterial culture.

Plasmids, as well as certain viruses, are extraordinarily useful tools for the molecular biologist, because they serve as gene-carrier molecules called cloning vectors. A basic procedure of recombinant DNA technology consists of joining a gene of interest to vector DNA to form a hybrid or recombinant molecule that is able to replicate in bacteria. Thus, cloning vectors contain genes for replication in bacteria. In addition, vectors generally carry antibiotic-resistance genes so that uninfected bacteria can be eliminated from the culture. In order to prepare a recombinant DNA molecule, a procedure is required for cutting cloning vectors and cellular DNA molecules at precise positions.

Figure 5. Introduction of a Plasmid Into Bacteria by Transformation



A nuclease is an enzyme that breaks the phosphodiester bonds that connect the nucleotide units in DNA or RNA. Restriction nucleases are powerful tools used in recombinant DNA technology because they cut DNA at specific sites. These enzymes are produced chiefly by bacterial species in which they serve to degrade invading foreign DNA within the bacterial cell. Most restriction enzymes recognize a specific sequence of four to six nucleotides in DNA and each will cut a long DNA double helix into a series of discrete pieces known as restriction fragments.

Table 1. Properties of Three Restriction Enzymes


Restriction Nuclease	Bacterial Source	Recognition Sequence and Mode of Cleavage
Hpa I	<i>Haemophilus parainfluenzae</i>	-G-T-T-A-A-C- -C-A-A-T-T-G- Cuts on lines of symmetry -G-T-T A-A-C- -C-A-A T-T-G- Blunt ends
EcoRI	<i>Escherichia coli</i>	-G-A-A-T-T-C- -C-T-T-A-A-G- Staggered Cleavage -G A-A-T-T-C -C-T-T-A-A G
Bam HI	<i>Bacillus amyloliquefaciens</i>	-G-G-A-T-C-C- -C-C-T-A-G-G- Staggered Cleavage -G G-A-T-C-C -C-C-T-A-G G Cohesive (Sticky) Ends

Typically, the restriction sites for a given enzyme are hundreds to thousands of base-pairs apart so that the fragments generated are hundreds to thousands of base-pairs long. More than 300 different restriction nucleases are now commercially available. General properties of three of these enzymes are given in Table 1. It should be noted that some restriction nucleases (e.g., EcoRI and Bam HI) produce a staggered cleavage that creates sticky, or cohesive, single-stranded ends on the cut molecules. These cohesive ends are very important in recombinant DNA procedures because they enable any two DNA fragments to be linked together by complementary base pairing at their ends, provided that they were generated with the same restriction enzyme.

Figure 6 illustrates one basic procedure by which cellular DNA from essentially any source can be amplified by recombinant DNA techniques. First, a plasmid is cleaved at a single site by a restriction nuclease, such as EcoRI, that produces cohesive ends on the plasmid DNA. The cellular DNA to be cloned is cleaved with the same enzyme, and fragments of the cellular DNA are annealed to the plasmid DNA by complementary base-pairing at the cohesive ends of the DNA molecules. The newly formed joints are sealed with an enzyme called DNA ligase, which forms covalent bonds between the ends of each DNA molecule. The recombinant DNA molecules are introduced into *E. coli* by transformation, and the bacteria are grown in the presence of an antibiotic. The hybrid plasmid can replicate in the dividing bacterial cells to produce an enormous number of copies of the original DNA fragment. At the end of the proliferation period, the hybrid plasmid molecules are purified from the bacteria. Copies of the original DNA fragments can then be recovered by cleavage of the recombinant plasmid with EcoRI.

Digestion of DNA from a vertebrate cell with EcoRI generates about 10^6 different DNA fragments. Thus, the DNA cloning procedure described above and outlined in Figure 6 gives rise to a large number of plasmids, each descended from a single hybrid DNA molecule. The most difficult step in the procedure is to identify the hybrid plasmid in this "library" that contains the inserted cellular DNA of interest. When a particular mRNA can be purified from a tissue, such as the mRNAs for the polypeptide chains of hemoglobin, the mRNA or a DNA copy of it can be used to identify its corresponding gene sequence in a recombinant library. In this procedure, the mRNA, or its DNA copy, is first labeled with a radioactive isotope. Under the appropriate conditions, the radioactive probe will preferentially stick or hybridize to the DNA clone of interest because of complementary base-pairing. By this procedure, the genes for many different proteins in DNA libraries have been identified.

Figure 6. Production of a Recombinant DNA Molecule



II. General Description of Agarose Gel Electrophoresis

Electrophoresis is the movement of charged particles in solution under the influence of an electric field. In the most common form of electrophoresis, the sample is applied to a stabilizing medium which serves as a matrix for the buffer in which the sample molecules travel. The agarose gel is a common type of stabilizing medium used for the electrophoretic separation of nucleic acids. A diagram of the essential components of an agarose electrophoretic system is shown in Figure 7. The agarose gel, containing preformed sample wells, is submerged in buffer within the electrophoretic gel cell. Samples to be separated are then loaded into the sample wells. Current from the power supply travels to the negative electrode (cathode), supplying electrons to the conductive buffer solution, gel and positive electrode (anode), thus completing the circuit.

At neutral pH, a molecule of DNA or RNA is negatively charged because of the negative charges on the phosphate backbone. Under these conditions, nucleic acids applied to sample wells at the negative electrode end of the gel migrate within pores of the gel matrix towards the positive electrode. The agarose gel serves as a molecular sieve in that its structure is similar to that of a sponge. The size of the pores in the gel are generally on the same order as the size of the DNA molecules that are being separated. As a result, large molecules move more slowly through the gel than smaller molecules. Thus, the method sorts the molecules according to size, since it relies on the ability of uniformly charged nucleic acids to fit through the pores of the agarose gel matrix.


III. Practical Aspects of Electrophoresis

Electrophoresis Equipment

The PROCELL Horizontal Electrophoresis unit is composed of an acrylic cell with central platform, platinum electrodes, four removable gel-casting trays, four sample well-forming combs and a safety lid with power cords. The four gels are made in the casting trays and then placed on the central platform of the electrophoresis cell. Each gel contains 8 separate sample wells. The experiments described below were designed such that each student uses four sample wells per experiment. Therefore, the experiments of 8 students can be analyzed in one electrophoretic run. If the students work in pairs, the system can be used by 16 students.

The Model MB-170 power supply is a general purpose electrophoresis power source. The unit produces a constant voltage output of 85 or 170 volts. Voltage selection is controlled by the switch located in the center of the front panel. The ammeter, also located on the front panel, permits the current to be monitored during an electrophoretic run. The unit can reach a maximum of 500 mAmp.

Figure 7. Components of a Horizontal Electrophoresis System.



Electrophoresis Chemicals

The Agarose Gel

Because the agarose gel is an ideal solid support for the separation of nucleic acids on the basis of size, it is used extensively for this purpose in the molecular biology laboratory. Agarose is a natural polysaccharide of galactose and 3,6-anhydrogalactose derived from agar, which, in turn, is obtained from certain marine red algae. Agarose gels are made by dissolving the dry polymer in boiling buffer, pouring the gels into casting trays and allowing them to set by cooling at room temperature. The resolving power of an agarose gel depends on the pore size, which

is dictated by the concentration of dissolved agarose. High percentage agarose gels (e.g., 3%) are used for the separation of small DNA molecules (10^2 - 10^3 base-pairs in length), while low percentage gels (e.g., 0.6%) are used for large molecules (10^5 - 10^6 base-pairs). The 1.2% and 1.5% agarose gels used in the experiments described in this manual are suitable for separating DNA molecules that range in length from about 200 to 10,000 base-pairs.

Electrophoresis Buffer

The buffer (0.04M Tris-Acetate-EDTA, pH 8.0) used in the experiments described below is a common buffer employed in the research laboratory for separating double-stranded DNA molecules by agarose gel electrophoresis. The buffer is provided as a 100-fold stock solution and should be diluted with distilled or deionized water prior to use.

Staining

Nucleic acids are not colored, and therefore it is necessary to make them visible in some way in order to determine their position in the agarose gel after electrophoresis. The most common method involves staining the nucleic acids with ethidium bromide, and then detecting the DNA or RNA-bound dye with ultraviolet lamps. However, ethidium bromide is a powerful mutagen (and probably a carcinogen) and UV light can cause serious eye and skin burns. Therefore, in the studies described below, the agarose gels are stained with the harmless dye, methylene blue.

The Sample Buffer

The DNA samples in the exercises described below are loaded into the wells of the agarose gel as 10-20% glycerol solutions. The viscous glycerol ensures that the samples will layer smoothly at the bottom of the sample wells. The sample buffer also contains the tracking dye bromophenol blue. As will be described below, this dye enables the investigator to follow the progress of an electrophoretic run.


Accessories

The accessories listed below are required to perform the experiments described in this manual. They are used for sample handling, and for the preparation and analysis of the agarose gels.

- *Glass test tubes (25ml)
- *Micro tubes (0.5ml)
- *Tube holders for the 0.5ml tubes
- *Tape
- *Gloves
- *Gel staining trays with lids
- *Macropipets (pipet-syringe)
- *Micropipets and micropipets

The micropipetors are an important accessory used for electrophoresis procedures. The micropipetor consists of a stainless steel plunger and 50 micropipets with calibration lines of 5 μ l*, 10 μ l, 15 μ l, 20 μ l, and 25 μ l as shown below.

A Diagram of the Micropipet Apparatus



*One milliliter (ml) = 1,000 microliters (μ l)

TO OPERATE THE MICROPIPETOR:

1. Insert the metal plunger into the end of the glass pipet that is opposite the calibration lines. The glass pipet can be held between your thumb and middle finger and the plunger operated with your index finger on the same hand.
2. Gently push down on the plunger until the plunger handle comes to rest on the pipet.
3. Hold the micropipet in a vertical position and place the filling end into the sample solution.
4. Draw the sample into the pipet to the appropriate calibration line by lifting up on the handle of the plunger assembly.
5. Carefully wipe excess liquid from the outer pipet surface with an absorbent tissue.
6. Direct the filling end of the micropipet into a tube and slowly eject the sample.
7. Rinse the pipet between samples by drawing up and expelling water three times from the pipetor.

Students should practice using these pipets prior to beginning the experiments.


IV. PROCEDURES FOR THE PREPARATION, ELECTROPHORESIS AND STAINING OF AGAROSE GELS

The exercises described in this manual were designed such that the samples of two students are analyzed on one agarose gel. If students work in pairs, four students will share one gel. Four agarose gels are electrophoresed simultaneously using the Procell and MB-170 power supply.

Pouring the Agarose Gels

1. Place the casting tray on a level work surface and place a precleaned glass slide into the gel support deck.
2. Seal both ends of the gel support deck with tape. The tape must be firmly pressed against the edges of the deck to ensure a tight seal.

CASTING TRAY ASSEMBLY



3. * With the micropipetor (pipet-syringe), dispense 1.5ml of electrophoresis buffer into a 25ml glass test tube and add 0.18 grams of agarose. The agarose can be weighed out directly on an appropriate balance. If a balance is not available, 0.18 grams of agarose can be estimated by filling a 0.5ml micro tube with agarose until two-thirds full. Gently swirl the glass tube until the agarose forms a suspension.

4. Place the test tube into a boiling water bath and allow the agarose suspension to come to a vigorous boil. After boiling for about two minutes, remove the test tube from the bath, stir gently with a glass rod, and cool at room temperature for about 2-3 minutes. At this time, the agarose solution should be absolutely clear.

*The melted agarose for the four gels (1.5ml per gel) can also be prepared in one operation by boiling 0.9g of agarose in 75ml of buffer in a 250-500ml flask over a bunsen burner or in a microwave oven. The flask should be rotated periodically during the heating process in order to prevent damage to the agarose.

5. Pour the melted agarose directly from the test tube onto the casting deck and return the test tube to the hot (but not boiling) water bath. The small amount of melted agarose left in the test tube will be used for sample application (see below). Insert the comb into the casting tray slots and push down gently on the top of the comb until resistance is encountered. The teeth of the comb will come to rest in the melted agarose about 0.2mm above the surface of the glass plate.
6. After the gel has cooled for at least 15 minutes, remove the tape strips and carefully lift the comb straight up and away from the casting tray. The gel is now ready for sample application. Gels can also be stored for up to one week before use. For gel storage, the comb is left in place and the tray containing the gel and comb is wrapped in plastic wrap and placed in the refrigerator.

Sample Application

Prior to sample application, place the following items on the laboratory bench in front of you.

- *Agarose gel
 - *Sample for electrophoresis
 - *Micropipetor and micropipets
 - *Absorbent tissue (e.g. Kleenex or Kimwipes)
 - *Small beaker of distilled or deionized water
 - *Melted agarose - Transfer the melted agarose from the large test tube in the hot water bath to a small tube. The small tube should then be placed in a beaker of hot water to ensure that the agarose remains in a liquid state.
1. Hold the micropipetor in a vertical position and place the filling end of the micropipet into the sample solution.
 2. Draw the sample into the pipet to the 15 μ l calibration line by lifting up on the handle of the plunger assembly.
 3. Wipe excess liquid from the outer pipet surface with an absorbent tissue.
 4. Carefully direct the filling end of the micropipet into the top of the sample well and slowly eject the 15 μ l of the sample well.
 5. Draw melted agarose into the micropipet to the 20 μ l calibration line, direct the filling end into the sample well, and slowly eject the agarose onto the sample until the well is full. Between 10-20 μ l of agarose are required to fill the well. The agarose will seal the sample in the sample well.
 6. Rinse the pipet by drawing up and expelling water three times from the pipetor.
 7. Wipe excess liquid from the outer pipet surface with an absorbent tissue.
 8. Repeat steps 1-7 to load each additional sample.

Electrophoresis

1. Transfer the four casting trays with gels to the central platform of the electrophoresis cell and position them such that the sample wells are closest to the black (negative) electrode. Upon electrophoresis, DNA will then migrate from the negative (black) towards the positive (red) electrode.
2. Place the gel tray stabilizing bar parallel to the long axis of the electrophoresis cell between the gel trays.
3. Slowly fill the electrophoresis chamber with electrophoresis buffer until the four gels are covered with a 1/4cm layer of buffer. Approximately 2.5 liters of buffer are required.
4. Place the electrophoresis cell lid in position.
5. With the power supply off, connect the cables from the cell to the power supply, red to red (positive) and black to black (negative).
6. Push the rocker switches on the power supply to "on" and "170V". The voltage will now remain constant at 170 volts during the run.
7. Unless otherwise indicated, electrophoresis until the bromophenol blue in the sample solution has migrated to within 1/4cm of the positive electrode end of the gel. At 170V, this takes approximately 50 minutes.
8. At the termination of the electrophoretic separation, shut off the power supply, disconnect the cables and remove the gel casting trays containing the gels.
9. The buffer should be emptied from the electrophoresis cell and stored under refrigeration in a separate container until the next electrophoretic run. The same buffer should be used for at least 3 electrophoretic separations. However, fresh buffer should be employed for the preparation of all agarose gels. The electrophoresis cell should be rinsed with deionized or distilled water and stored in an upright position.

Staining

Normal Procedure (Overnight Staining)

1. Carefully slide the agarose gel out of the casting tray and off of the slide and place the gel in a staining dish. (Note: The glass slide should not be placed in the staining dish.)
2. Dilute the gel stain concentrate 1000-fold with distilled water.
3. Cover the gel with about 100ml of staining solution, making certain that the agarose does not stick to the dish and all gel surfaces are exposed to the stain.

4. Place the staining tray in the refrigerator and allow 3-18 hours for staining.
5. Decant and discard the stain, rinse the gel and dish with distilled water and add about 100ml of distilled water.
6. Change the water after about 10-20 minutes and hold the staining dish over a light source such as a desk lamp or light box. Note the position of the dark blue DNA bands against the light blue gel background. The gel can be stored in water in the refrigerator for a few days. After about a week or two in the refrigerator, the dye will diffuse out of the gel and the DNA bands will no longer be visible.

Rapid Procedure

1. Dilute the gel stain concentrate 1000-fold and stain the gels for 30-60 minutes in 100mls at 37°C.
2. Decant and discard the stain, rinse the gel and dish with water and add about 100ml of distilled water.
3. Change the water after about 20 minutes and again after an additional 20 minutes. DNA bands can be seen during these destaining steps.
4. When the background stain has been reduced sufficiently, hold the staining dish over a light source and carefully observe the stained DNA bands in the gel.

Gel Storage

1. The gel can be stored in a sealed plastic bag (3x3") with a few ml's of distilled water for up to one month in the refrigerator.
2. For long term storage, place the gel on a glass slide and smooth with a gloved index finger to eliminate air bubbles between the gel and the slide. Allow the gel to dry onto the slide at room temperature for 3-4 days. Cover the dry gel film and glass slide with saran wrap.

V. Suggested Reading and Reference for Part A

Lewin, B. *Genes V*. Oxford University Press, Oxford New York Tokyo, 1994.

Jones, P. *Gel Electrophoresis: Nucleic Acids*. Chichester, West Sussex, UK, New York: Wiley 1995.